

Comprehensive Helicase Solutions: Protein Production, Biochemical, Biophysical Assays, Selectivity Panel & Safety Assessment

Jean Bernatchez¹, Romain Merceron³, Thibault Anna⁴, Vanessa Jahnke⁴, Vanessa Porkolab⁴, Sophie Raynal⁴, Emilie Desfosses⁴, Samuel Ramirez¹, Fabien Tillier¹, Celine Legros²

1. Eurofins Discovery | USA, 2. Eurofins Discovery | France, 3. Eurofins DiscoverX Products France | France, 4. Eurofins Cerep | France


Introduction

Helicases are specialized enzymes, essential for unwinding DNA-DNA or RNA-DNA duplexes and required for genome maintenance and cellular homeostasis. Their mutation and/or overexpression are closely associated with many disorders including tumor development and progression. These cancer-specific vulnerabilities can thus be exploited to selectively target tumors. Eurofins Discovery's experts have worked together to build integrated and comprehensive Helicase Drug Discovery solutions from protein production to safety profiling, to quickly and efficiently identify best drug candidates.

The first part of the poster is dedicated to high-quality helicase proteins compatible with biochemical and biophysics applications, allowing profiling studies, high-throughput screening (HTS) and fragment-based screening (FBS). The second part of the poster shows how to use the different products and services to study and de-risk a helicase inhibitor.


Protein Production & QC

Eurofins DiscoverX® generates high-quality protein, as exemplified by the production of truncated PolQ (Polymerase Theta), DNA helicase RECQ1 and RNA helicase RIG-I, exhibiting a purity exceeding 90%.



Biophysics Characterization

These proteins are biophysics SPR compatible. We were able to measure affinity binding to ADP with a K_D of 19.5 μM for truncated PolQ (helicase domain), 23.5 μM for RECQ1 DNA helicase and 19 μM for RIG-I RNA helicase.



Biochemical Enzyme Activity

Additionally, ATPase and unwinding activity assays were developed to demonstrate that proteins are enzymatically functional. The example shows the principle of ATPase activity assay and a linearity of signal between ATPase activity and enzyme concentration.





Figure 3. ATPase activity. **A.** The ATPase activity was performed using DNA and RNA double strand with ADP Hunter Plus assay kit (#90-0083, Eurofins DiscoverX). Activity was monitored by fluorescence through the accumulation of ADP. **B.** PolQ (helicase domain) protein activity – enzyme titration. **C.** RECQ1 protein activity – enzyme titration. **D.** RIG-I protein activity – enzyme titration.

Biochemical Characterization of a PolQ Inhibitor

To validate our materials and methods, we used the specific inhibitor PolQi2. We first determined the characteristics of the enzyme (K_M substrate and K_M ATP) before measuring IC₅₀ of PolQi2 (4nM), in line with the literature. We also characterized the mode of inhibition, and demonstrated that PolQi2 is an allosteric inhibitor.



PolQ Inhibitor Selectivity Panel


To check the selectivity of PolQi2 against other helicases from the same super family, we used our selectivity panel and confirmed that PolQi2 is specific to the PolQ enzyme.

	RECQ1	BLM/RECQ2WRN/RECQ3	RTS/RECQ4	RECQ5	FANCJ	PolQ	DHX58/LGP2	RIG-1/DDX58	MDA5/IFIH1
PolQi2						4nM			
ML216		2μM	3μM	8μM	1.5μM	4μM	3.5μM		5μM
Tyrphostin AG538	1μM	2μM	1.5μM	2.5μM	1.5μM	3.5μM	4μM	5μM	1μM

Table 1. PolQi2 a selective inhibitor of PolQ. Selectivity testing of 3 helicases inhibitors, PolQi2, ML216 and Tyrphostin. AG538, in the helicases ATPases activity panel. PolQi2 is highly selective for PolQ in comparison to poorly selective inhibitors

SafetyScreen 18 Core Panel

According to Brennan et al. (2024), an early de-risking is recommended to select or eliminate chemical series with adverse safety alerts. We therefore performed our SafetyScreen 18 Core Panel on PolQi2 and found no red alerts using our hazard metric tool.



Helicase Assays - Overview

Eurofins DiscoverX					Eurofins Discovery Services			
Helicase	Substrate	Tags	Purity	Protein Item Number	ATPase Item Number	Unwinding Item Number	SPR	MST TRIC / Spectral Shift
RECQ1	DNA	6HisAvi	FL - >95%	18-005	5763	Ongoing	Validated	Ongoing
BLM/RECQ2	DNA	MBP, 6HisAvi	FL - 92%	18-010	5751	5752	n/a	n/a
WRN/RECQ3	DNA	MBP, 6HisAvi	FL - 88%	18-009	5749	5750	n/a	n/a
RTS/RECQ4	DNA	MBP, 6HisAvi	FL - 90%	18-006	5765	Ongoing	n/a	n/a
RECQ5	DNA	6HisAvi	FL - 90%	18-007	5753	Ongoing	Validated	Ongoing
FANCJ	DNA	MBP, 6HisAvi	FL - 85%	18-001	5755	Ongoing	n/a	n/a
PolQ (1-894)	DNA	6HisAvi	Helicase - >95%	18-004	5761	Ongoing	Validated	Validated
DHX58/LGP2	RNA	6HisAvi	FL - 74%	18-002	5757	Ongoing	Validated	Ongoing
MDA5/IFIH1	RNA	6HisAvi	FL - >95%	18-003	5759	Ongoing	Validated	Validated
RIG-I/DDX58	RNA	6HisAvi	FL - >95%	18-008	5767	Ongoing	Validated	Validated

Table 2. Current products and assays available. Nine full length and one truncated recombinant helicase proteins are available. Off-the-shelf biochemistry and biophysics assays currently in our helicase offer.

Summary

Relying on the expertise of the Eurofins DiscoverX team to produce high-quality proteins and ADP measurement kit, Eurofins Discovery has successfully developed a set of 10 DNA & RNA helicases (and soon a set of 8 new RNA helicases), illustrated here with the example of Polymerase Theta, PolQ. The access of a selectivity panel opens the door to new drug discovery challenges and the metric tool associated with our Safety Core panel 18 offer a valuable support to early safety assessment.

The combination of the Eurofins Discovery helicase platform with other pre-clinical drug discovery capabilities, supported by our experienced team of Scientific Experts and Project Managers, gives you the best chance to quickly and efficiently identify best drug candidates.